

FOLEY CATHETER

UROLOGY

- Symmetrical large capacity balloon ensuring a straight tip for proper flow and good sphincter action that prevent bladder leakage.
- Manufactured from natural rubber coated with silicone.
- Coned distal end provided with burr free eyes for atraumatic intubation.
- Hard valve for easy inflation and deflation of balloon.


▪ Optimally sized drainage eye ensure effective drainage with enhance tip strength.

▪ Strong pressure-resistant balloon for effective anchoring at bladder neck.

▪ Thick and rigid funnel ensure secure hold with drainage bag's connector.




▪ Self-sealing inflation hard valve ensure balloon is inflated when in use.

▪ Strong catheter shaft for effective insertion and removal yet flexible to maximize patient's comfort.



Standard Foley Catheter Silicone Coated, Sterile

Various tips configuration for different catheterization needs

	STANDARD
	COUVELAIRE
	DUFOUR


All Silicone Foley Catheter, Sterile




Standard Foley Catheter Silicone Coated, Sterile




Premium Foley Catheter Gold, Siliconised, Sterile



Nelaton Premium Foley Catheter Gold, Siliconised, Sterile



Premium Hematuria Reinforced Foley Catheter Gold, Siliconised, Sterile



Ideal Healthcare Sdn. Bhd.(484861-V)
No. 70 & 71, Jalan Sungai Tukang 2/1, Kawasan Perusahaan Sungai Tukang,
08000 Sungai Petani, Kedah, Malaysia.
Tel : 6 (04) - 422 7775 / 422 8775 Fax : 6 (04) - 422 6775
Email : info@idealcare.com.my

E. L. NO.: MDA-0003-KD314

COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 13485 =

www.idealcare.com.my



idealcareMY

IHSB-CAT-BDD-TDR-URO-02

All Silicone Foley Catheter, Sterile

2-Way Paediatric (HV)		2-Way Standard (HV)		3-Way Standard (HV)		Size (Ch/Fr)
1-5ml/cc	3-5ml/cc	5-15ml/cc	30-45ml/cc	5-15ml/cc	30-45ml/cc	
IC43150603	-	-	-	-	-	06
-	IC43150803	-	-	-	-	08
-	IC43151003	-	-	-	-	10
-	-	IC41151205	-	-	-	12
-	-	IC41151405	-	-	-	14
-	-	IC41151605	IC41151630	-	IC45151630	16
-	-	IC41151805	IC41151830	IC45151805	IC45151830	18
-	-	IC41152005	IC41152030	IC45152005	IC45152030	20
-	-	IC41152205	IC41152230	IC45152205	IC45152230	22
-	-	IC41152405	IC41152430	IC45152405	IC45152430	24
-	-	IC41152605	IC41152630	IC45152605	IC45152630	26

Note: HV = Hard Valve

Standard Foley Catheter Silicone Coated, Sterile

2-Way Paediatric (HV)		2-Way Standard (HV)		3-Way Standard (HV)		Size (Ch/Fr)
3-5ml/cc	5-15ml/cc	20-30ml/cc	30-45ml/cc	5-15ml/cc	30-45ml/cc	
IC23120605	-	-	-	-	-	06
IC23120805	-	-	-	-	-	08
IC23121005	-	-	-	-	-	10
-	IC21111210	IC21111230	-	-	-	12
-	IC21111410	IC21111430	-	-	-	14
-	IC21111610	-	IC21111630	IC25121610	IC25121630	16
-	IC21111810	-	IC21111830	IC25121810	IC25121830	18
-	IC21112010	-	IC21112030	IC25122010	IC25122030	20
-	IC21112210	-	IC21112230	IC25122210	IC25122230	22
-	IC21112410	-	IC21112430	IC25122410	IC25122430	24
-	IC21112610	-	IC21112630	IC25122610	IC25122630	26
-	IC21112810	-	IC21112830	-	-	28
-	IC21113010	-	IC21113030	-	-	30

Note: HV = Hard Valve

- ✓ Disposable.
- ✓ Sterile.
- ✓ Easy to use.
- ✓ Individually packed.
- ✓ Packaging: 10pcs/box, 50 boxes / carton.

Premium Foley Catheter Gold, Siliconised, Sterile

2-Way Premium Paediatric Foley Catheter	2-Way Premium Foley Catheter		3-Way Premium Foley Catheter		Size (Ch/Fr)
3-5ml/cc	5-15ml/cc	30-45ml/cc	5-15ml/cc	30-45ml/cc	
GD23120605	-	-	-	-	06
GD23120805	-	-	-	-	08
GD23121005	-	-	-	-	10
-	GD21111210	GD21111230	-	-	12
-	GD21111410	GD21111430	-	-	14
-	GD21111610	GD21111630	GD25121610	GD25121630	16
-	GD21111810	GD21111830	GD25121810	GD25121830	18
-	GD21112010	GD21112030	GD25122010	GD25122030	20
-	GD21112210	GD21112230	GD25122210	GD25122230	22
-	GD21112410	GD21112430	GD25122410	GD25122430	24
-	GD21112610	GD21112630	GD25122610	GD25122630	26

2-Way Premium Hematuria Reinforced Foley Catheter Gold, Siliconised, Sterile

Couvelaire		Dufour		Size (Ch/Fr)
30-50 ml/cc	60-80 ml/cc	30-50 ml/cc	60-80 ml/cc	
GD21421830	-	-	-	18
GD21422030	GD21422060	GD21522030	GD21522060	20
GD21422230	GD21422260	GD21522230	GD21522260	22
GD21422430	GD21422460	GD21522430	GD21522460	24
GD21422630	GD21422660	GD21522630	GD21522660	26

3-Way Premium Hematuria Reinforced Foley Catheter Gold, Siliconised, Sterile

Couvelaire		Dufour		Size (Ch/Fr)
30-50 ml/cc	60-80 ml/cc	30-50 ml/cc	60-80 ml/cc	
GD25421830	-	-	-	18
GD25422030	GD25422060	GD25522030	GD25522060	20
GD25422230	GD25422260	GD25522230	GD25522260	22
GD25422430	GD25422460	GD25522430	GD25522460	24
GD25422630	GD25422660	GD25522630	GD25522660	26

Nelaton Premium Foley Catheter Gold, Siliconised, Sterile

Catalogue No.	Size
GD36101200	12
GD36101400	14
GD36101600	16
GD36101800	18
GD36102000	20
GD36102200	22
GD36102400	24
GD36102600	26