



**IDEALCARE<sup>®</sup>**

*.....the quality you can trust*

PRODUCT  
CATALOGUE



[www.idealcare.com.my](http://www.idealcare.com.my)

## **VISION**

To be a one-stop centre for hospitals, healthcare centers, clinics and pharmacies. Already a market leader in Malaysia and with the continuous improvement attitude it will become a market leader regionally as well as internationally.

## **MISSION**

Support and enhance the market leadership of the Disposable Product Businesses through the distribution of the highest quality products at superior service levels and most competitive prices.

The foundation for this effort will be the full partnership of our employee in the attainment of our goals through open communications and total involvement in our business.

## **QUALITY POLICY**

"Ideal Healthcare Group of Companies are committed to providing the highest standard quality products and services that will meet or exceed our customers' expectation, comply with regulatory requirements, and to maintain the effectiveness of the quality management system."



# CORPORATE PROFILE

Ideal Healthcare Sdn Bhd (IHSB) formally known as Mutiara Medical was established in 1999. Founded and steered by Dato Dr Haji Haminnuddin Bin Haji Abd. Hamid as the Chief Executive, IHSB specializes its corporate and operational activities in the scope of manufacturing, assemblies and distribution of medical disposable products particularly in the field of Urology, Anesthesia, Respiratory, Infusion / Transfusion Therapy and Procedure Sets under the registered trademark of IDEALCARE®

At IHSB, we focus on serving the need of our customers by offering superior quality products and services at extremely competitive price. We are also committed in ensuring safety and performance of our product through effective implementation Quality Management System in conformance to ISO 13485:2016, Good Distribution Practice for Medical Device (GDPMD) and compliance against the National/International Medical Devices Regulation.

Over the years, IHSB consistently improve the organization value and market presence by successfully establishing itself as noticeable competitor to many Multi-National Corporation in the medical device industry both domestically and globally. As an acknowledgement against our reputation in quality and customer oriented service, IHSB have regularly been selected as substantial supplier to the Ministry of Health Malaysia since 2005 to supply the entire spectrum of its products to the government hospitals, healthcare centers, clinics and institutions throughout Malaysia. IHSB have also penetrated the market of more than 20 countries in the South East Asia, Middle East and South Pacific to further evident the global market recognition and acceptance towards our brands, product and services.

Moving forward, IHSB have recently completed the development of its latest state of the art facility to accommodate the growing worldwide demand on our products. Well-known by the name of IHSB Centre of Excellence (COE), these facility is entirely run by automation and equipped with R&D laboratory to demonstrate our persistent commitment and vision to progress beyond expectation for the future need of our consumers. Today, with more than 20 years of industry experience, IHSB continue to foster the tradition of excellence product and service performance to ensure IDEALCARE® will become the world renown brand of choice in years to come.



**DATO' DR HJ HAMINNUDDIN BIN HJ ABD HAMID**  
Group Founder & Chief Executive

# Contents

## Anaesthesia

Endotracheal Tube	1-2
Tracheostomy Tube	3
Laryngeal Airway Mask	4
Suction Catheter	4
Nasopharyngeal Airway	4

## Respiratory Care

Resuscitator	5
Oxygen Mask with Tubing	5
Nasal Oxygen Cannula	5
Oxygen Mask High Concen. Non Rebreathing Mask	5
Nebulizer Mask with Chamber & Tubing	6
Oxygen Air Entrainment	6
Heat Moisture Exchange (HME) Filter	6

## Urology

Foley Catheter	7-8
Adult Urine Collecting Bag	9
Paediatric Urine Collector	9
Male External Catheter	9
Urine Collecting Bag with Measured Volume Meter	10
Urine Bag Hanger	10

## Infusion & Transfusion

Intravenous Catheter with Injection Valve, Quick Flashback	11
Intravenous Catheter with Injection Port & Wings	12
Intravenous Catheter with Injection Valve	13
Intravenous Catheter Without Port & Without Wings	14
Blood Giving Set with Airvent (Blood Transfusion Set)	15
Intravenous Administration Set	15-16
Microdrip Precision Solution Administration Set	17
Peritoneal Dialysis Transfusion Set	17
Peritoneal Dialysis Catheter Set	17
Extension Tube With Male Female Luer Lock	18
CVP Manometer	18
Extension Tube with 'T' Connector	18

## Gastroenterology

Ryle's Tube	19
Infant Feeding Tube	19
Umbilical Catheter	19
Mucus Extractor	19

## Surgery & Wound Drainage

Under Water Seal Pleural Drainage System	20
Connecting Tube with Clamp & 'Y' Connector	20
Perforated Redon Drain Catheter	20
Close Wound Suction Unit	21
Suction Connecting Tube	21
Yankaur Suction Set	21
Yankaur Suction Handle (Standard tip & Crown tip)	21
Colostomy Bag	22
Thoracic Trocar Catheter	22
Chest Drainage Catheter	22

## Pathology

Petri Dishes	23
Proctoscope	23
Urine Specimen Container	23
Sputum Cup	23

## Obstetrics & Gynaecology

Amniotic Membrane Perforator	24
Vaginal Speculum	24
Cervical Brush	24
Gynecological Brush with Tip Ball End	24

## Procedure & Clinical Sets

Basic Dressing Set (Premium)	24
Basic Dressing Set (Standard)	25
Basic Dressing Set (Haemo)	25
Vaginal Examination Set	25
Catheterization Set	25
Delivery Set	26
Injection Set	26
Pap Smear Test Kit	26
Toilet & Suture Set	26
Top Tail/Baby Cleaning Set	26
Eye Dressing Set	27
Cord Dressing Set	27
Suction Set	27
Blood Taking Set	27
Swabbing Set	27

## Procedure Set Components

Kidney Dish	28
Forceps	28
Umbilical Cord Clamp	28

## Wound Care & Dressing

Cotton Wool Roll	28
Combine Dressing Roll	28
White Elastic Adhesive Bandage	28
Gauze Ribbon	29
Lint Absorbent	29
Plaster Extension	29
Skin Color Crepe Bandage	29
Skin Traction Kit	29
Cotton Wool Balls	29
Transparent Film Dressing	30
Transparent Infusion Set Dressing	30
Eye Pad	30
Paraffin Gauze Dressing	30
Soft Paper Surgical Tape	30
Surgical Drape	30
Alcohol Swab	31
Wooden Applicator (Orange Stick)	31
Wooden Tongue Depressor	31


## Rehabilitation

Wheelchair	31
------------	----




## Endotracheal Tube

- Safe & effective airway management.
- Suitable for general & emergency applications.
- User friendly for anaesthetist & patients.
- Easy insertion and less trauma.



- Split-resistant radio-opaque line ensure tube integrity when cut to required length and connector reinserted. Quality material results in **clear radiopacity**, thus facilitating verification of intubated position.
- Smooth, flexible and kink resistant inflation line ensures for **effective inflation & deflation process**.
- Non-return valve to ensure continual integrity and size marking on pilot balloon for better **product identification**.
- Double-ring marks help guide intubation through laryngoscopy.
- Soft bull-nose tip enables smooth **atraumatic intubation**.
- High quality cuff for **effective sealing**.

Wire Reinforced Endotracheal Tubes Cuffed With Stylet



Endotracheal Tube Regular Plain Uncuff



### Common Features













Murphy eye at distal end for atraumatic intubation and extubation.



Radiopaque line on the tube to facilitate easy location of the tube.



Proximal end fitted with standard 15mm connector.


	CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
 REGULAR HVLP CUFFED	IC 61110030	3.0 mm	IC 61110070	7.0 mm
	IC 61110035	3.5 mm	IC 61110075	7.5 mm
	IC 61110040	4.0 mm	IC 61110080	8.0 mm
	IC 61110045	4.5 mm	IC 61110085	8.5 mm
	IC 61110050	5.0 mm	IC 61110090	9.0 mm
	IC 61110055	5.5 mm	IC 61110095	9.5 mm
	IC 61110060	6.0 mm	IC 61110100	10.0 mm
	IC 61110065	6.5 mm	-	-
 REGULAR SVLP CUFFED	IC 61150030	3.0 mm	IC 61150065	6.5 mm
	IC 61150035	3.0 mm	IC 61150070	7.0 mm
	IC 61150040	4.0 mm	IC 61150075	7.5 mm
	IC 61150045	4.5 mm	IC 61150080	8.0 mm
	IC 61150050	5.0 mm	IC 61150085	8.5 mm
	IC 61150055	5.5 mm	IC 61150090	9.0 mm
	IC 61150060	6.0 mm	-	-
	IC 61214030	3.0 mm	IC 61214065	6.5 mm
 WIRE REINFORCED HVLP CUFFED	IC 61214035	3.5 mm	IC 61214070	7.0 mm
	IC 61214040	4.0 mm	IC 61214075	7.5 mm
	IC 61214045	4.5 mm	IC 61214080	8.0 mm
	IC 61214050	5.0 mm	IC 61214085	8.5 mm
	IC 61214055	5.5 mm	IC 61214090	9.0 mm
	IC 61214060	6.0 mm	IC 61214095	9.5 mm
	IC 61214130	3.0 mm	IC 61214165	6.5 mm
	IC 61214135	3.5 mm	IC 61214170	7.0 mm
 WIRE REINFORCED HVLP CUFFED WITH STYLET	IC 61214140	4.0 mm	IC 61214175	7.5 mm
	IC 61214145	4.5 mm	IC 61214180	8.0 mm
	IC 61214150	5.0 mm	IC 61214185	8.5 mm
	IC 61214155	5.5 mm	IC 61214190	9.0 mm
	IC 61214160	6.0 mm	IC 61214195	9.5 mm
	IC 61213040	4.0 mm	IC 61213070	7.0 mm
	IC 61213045	4.5 mm	IC 61213075	7.5 mm
	IC 61213050	5.0 mm	IC 61213080	8.0 mm
 ORAL CURVED SVLP CUFFED	IC 61213055	5.5 mm	IC 61213085	8.5 mm
	IC 61213060	6.0 mm	IC 61213090	9.0 mm
	IC 61213065	6.5 mm	-	-
	IC 61313045	4.5 mm	IC 61313070	7.0 mm
	IC 61313050	5.0 mm	IC 61313075	7.5 mm
	IC 61313055	5.5 mm	IC 61313080	8.0 mm
	IC 61313060	6.0 mm	IC 61313085	8.5 mm
	IC 61313065	6.5 mm	IC 61313090	9.0 mm
 NASAL CURVED SVLP CUFFED	IC 61134025	2.5 mm	IC 61134065	6.5 mm
	IC 61134030	3.0 mm	IC 61134070	7.0 mm
	IC 61134035	3.5 mm	IC 61134075	7.5 mm
	IC 61134040	4.0 mm	IC 61134080	8.0 mm
	IC 61134045	4.5 mm	IC 61134085	8.5 mm
	IC 61134050	5.0 mm	IC 61134090	9.0 mm
	IC 61134055	5.5 mm	IC 61134095	9.5 mm
	IC 61134060	6.0 mm	-	-
 WIRE REINFORCED PLAIN UNCUFF	IC 61130020	2.0 mm	IC 61130065	6.5 mm
	IC 61130025	2.5 mm	IC 61130070	7.0 mm
	IC 61130030	3.0 mm	IC 61130075	7.5 mm
	IC 61130035	3.5 mm	IC 61130080	8.0 mm
	IC 61130040	4.0 mm	IC 61130085	8.5 mm
	IC 61130045	4.5 mm	IC 61130090	9.0 mm
	IC 61130050	5.0 mm	IC 61130095	9.5 mm
	IC 61130055	5.5 mm	IC 61130100	10.0 mm
	IC 61130060	6.0 mm	-	-
	IC 61233030	3.0 mm	IC 61233065	6.5 mm
 ORAL CURVED PLAIN UNCUFF	IC 61233035	3.5 mm	IC 61233070	7.0 mm
	IC 61233040	4.0 mm	IC 61233075	7.5 mm
	IC 61233045	4.5 mm	IC 61233080	8.0 mm
	IC 61233050	5.0 mm	IC 61233085	8.5 mm
	IC 61233055	5.5 mm	IC 61233090	9.0 mm
	IC 61233060	6.0 mm	-	-
	IC 61333030	3.0 mm	IC 61333065	6.5 mm
	IC 61333035	3.5 mm	IC 61333070	7.0 mm
 NASAL CURVED PLAIN UNCUFF	IC 61333040	4.0 mm	IC 61333075	7.5 mm
	IC 61333045	4.5 mm	IC 61333080	8.0 mm
	IC 61333050	5.0 mm	IC 61333085	8.5 mm
	IC 61333055	5.5 mm	IC 61333090	9.0 mm
	IC 61333060	6.0 mm	-	-
	IC 64200006	6Fr (2.0 mm)	Length: 305 mm	-
	IC 64200010	10Fr (3.3 mm)	Length: 390 mm	-
	IC 64200014	14Fr (4.7 mm)	Length: 420 mm	-
 INTUBATION STYLET				



# Tracheostomy Tube

## Cuffed

## Uncuff



CATALOGUE NO.	SIZE
IC62520040	4.0 mm
IC62520045	4.5 mm
IC62520050	5.0 mm
IC62520055	5.5 mm
IC62520060	6.0 mm
IC62520065	6.5 mm
IC62520070	7.0 mm
IC62520075	7.5 mm
IC62520080	8.0 mm
IC62520085	8.5 mm
IC62520090	9.0 mm
IC62520095	9.5 mm
IC62520100	10.0 mm

CATALOGUE NO.	SIZE
IC62540030	3.0 mm
IC62540035	3.5 mm
IC62540040	4.0 mm
IC62540045	4.5 mm
IC62540050	5.0 mm
IC62540055	5.5 mm
IC62540060	6.0 mm
IC62540065	6.5 mm
IC62540070	7.0 mm
IC62540075	7.5 mm
IC62540080	8.0 mm
IC62540085	8.5 mm
IC62540090	9.0 mm
IC62540095	9.5 mm
IC62540100	10.0 mm

## Laryngeal Airway Mask

- Available in medical-grade PVC & 100% silicone.
- Safe & effective alternative to the endotracheal tube and face mask.
- Suitable for general & emergency applications.
- Reusable type can be resterilized.



### Disposable (PVC)

CATALOGUE NO.	SIZE	PATIENT WEIGHT
IC6310D	1.0 mm	Neonate < 5 kgs
IC6315D	1.5 mm	Infant, 5-10 kgs
IC6320D	2.0 mm	Infant, 10-20 kgs
IC6325D	2.5 mm	Large Child, 20-30 kgs
IC6330D	3.0 mm	Child, 30-50 kgs
IC6340D	4.0 mm	Adult, 50-70 kgs
IC6350D	5.0 mm	Large Adult, 70-100 kgs


### Reusable (Silicone)

CATALOGUE NO.	SIZE	PATIENT WEIGHT
IC6410R	1.0mm	Neonate < 5 kgs
IC6415R	1.5mm	Infant, 5-10 kgs
IC6420R	2.0mm	Infant, 10-20 kgs
IC6425R	2.5mm	Large Child, 20-30 Kgs
IC6430R	3.0mm	Child, 30-50 Kgs
IC6440R	4.0mm	Adult, 50-70 kgs
IC6450R	5.0mm	Large Adult, 70-100 kgs

## Suction Catheter

- Suitable for removal of secretion from mouth, trachea and bronchial tubes.
- Manufactured from non-toxic, non-irritant medical grade PVC.
- Frozen surface catheter tubing for smooth intubation.

CATALOGUE NO.	SIZE	COLOR CODE
IC2030-05	5Fr	Grey
IC2030-06	6Fr	Green
IC2030-08	8Fr	Blue
IC2030-10	10Fr	Black
IC2030-12	12Fr	White
IC2030-14	14Fr	Dark Green
IC2030-16	16Fr	Orange
IC2030-18	18Fr	Red



Suction Catheter with Finger Tip Control



## Nasopharyngeal Airway

- Suitable for suctioning and provide alternative way when oropharyngeal airway is occluded.
- Made of Medical Grade PVC readily conforms to the nasopharyngeal anatomy.
- Interior walls with channel grooves, facilitate easy insertion of catheter.

CATALOGUE NO.	SIZE
IC2034-6.5	6.5 mm
IC2034-7.0	7.0 mm
IC2034-7.5	7.5 mm
IC2034-8.0	8.0 mm
IC2034-8.5	8.5 mm



## Manual Resuscitator

- Available in medical-grade PVC & silicone.
- Autoclavable and made of medical grade silicone material.
- Easy to disassemble for cleaning or disinfecting.

### Silicone

CATALOGUE NO.	SIZE
IC2103	Adult
IC2203	Child
IC2303	Infant

### PVC

CATALOGUE NO.	SIZE
IC2501	Adult
IC2601	Child
IC2701	Infant



## Oxygen Mask with Tubing

- Designed for convenient oxygen therapy.
- Moulded face mask is made of medical grade PVC.

CATALOGUE NO.	SIZE
IC2020A	Elongated Adult
IC2020P	Elongated Paediatric
IC2020I	Elongated Infant
IC2020N	Elongated Neonate

## Nasal Oxygen Cannula

- Easy application and efficient administration of oxygen for patient's comfort.

### Non-sterile

CATALOGUE NO.	SIZE
IC2016A	Adult
IC2016P	Paediatric
IC2016I	Infant
IC2016N	Neonate

### Sterile

CATALOGUE NO.	SIZE
IC2017A	Adult
IC2017P	Paediatric



## Oxygen Mask High Concentration Non Rebreathing Mask

- Equipped with soft, transparent and odour free mask.
- Light weight reservoir improve the oxygen efficiency.

CATALOGUE NO.	SIZE
IC1110A	Adult
IC1112P	Paediatric

## Nebulizer Mask with Chamber & Tubing

- Soft and clear mask with anatomical form is ideal for long term use.
- Intergrated nose bridge for extra comfort.

CATALOGUE NO.	SIZE
IC1120	Elongated Adult
IC1122	Elongated Paediatric
IC1123	Standard Neonate



## Oxygen Air Entrainment

- Made of medical grade PVC.
- Adjustable nose clip.

Adjusted Oxygen Concentration	Color Code
24% (Blue)	
28% (Yellow)	
35% (Dark Green)	
40% (Pink)	
60% (Red)	



CATALOGUE NO.	SIZE
IC7411	Adult
IC7412	Paediatric

## Heat Moisture Exchange (HME) Filter

- Made of medical grade material.
- Mainly used for filtering bacteria.

CATALOGUE NO.	SIZE
IC158000	Standard






# Foley Catheter

## 2-Way Standard Foley Catheter Silicone Coated

- Symmetrical large capacity balloon ensuring a straight tip for proper flow and good sphincter action that prevent bladder leakage.
- Manufactured from natural rubber coated with silicone.
- Coned distal end provided with burr free eyes for atraumatic intubation.
- Hard valve for easy inflation and deflation of balloon.



- Thick and rigid funnel ensure secure hold with drainage bag's connector.
- Self-sealing inflation hard valve ensure balloon is inflated when in use.
- Strong catheter shaft for effective insertion and removal yet flexible to maximize patient's comfort.

- Optimally sized eyelet ensure effective drainage with enhance tip strength.
- Strong pressure-resistant balloon for effective anchoring at bladder neck.

### All Silicone Foley Catheter

### Premium Foley Catheter Gold, Siliconised



2-Way All Silicone Foley Catheter

2-Way Premium Foley Catheter Gold, Siliconised

2-Way Paediatric (HV)		2-Way Standard (HV)		3-Way Standard (HV)		Size (Ch/Fr)	All Silicone Foley Catheter
1-5ml/cc	3-5ml/cc	5-15ml/cc	30-45ml/cc	5-15ml/cc	30-45ml/cc		
IC43150603	-	-	-	-	-	06	Note: HV = Hard Valve
-	IC43150803	-	-	-	-	08	
-	IC43151003	-	-	-	-	10	
-	-	IC41151205	-	-	-	12	
-	-	IC41151405	-	-	-	14	
-	-	IC41151605	IC41151630	-	IC45151630	16	
-	-	IC41151805	IC41151830	IC45151805	IC45151830	18	
-	-	IC41152005	IC41152030	IC45152005	IC45152030	20	
-	-	IC41152205	IC41152230	IC45152205	IC45152230	22	
-	-	IC41152405	IC41152430	IC45152405	IC45152430	24	
-	-	IC41152605	IC41152630	IC45152605	IC45152630	26	

2-Way Paediatric (HV)		2-Way Standard (HV)		3-Way Standard (HV)		Size (Ch/Fr)	Standard Foley Catheter Silicone Coated
3-5ml/cc	5-15ml/cc	20-30ml/cc	30-45ml/cc	5-15ml/cc	30-45ml/cc		
IC23120605	-	-	-	-	-	06	Note: HV = Hard Valve
IC23120805	-	-	-	-	-	08	
IC23121005	-	-	-	-	-	10	
-	IC21111210	IC21111230	-	-	-	12	
-	IC21111410	IC21111430	-	-	-	14	
-	IC21111610	-	IC21111630	IC25121610	IC25121630	16	
-	IC21111810	-	IC21111830	IC25121810	IC25121830	18	
-	IC21112010	-	IC21112030	IC25122010	IC25122030	20	
-	IC21112210	-	IC21112230	IC25122210	IC25122230	22	
-	IC21112410	-	IC21112430	IC25122410	IC25122430	24	
-	IC21112610	-	IC21112630	IC25122610	IC25122630	26	
-	IC21112810	-	IC21112830	-	-	28	
-	IC21113010	-	IC21113030	-	-	30	

2-Way Premium Paediatric Foley Catheter		2-Way Premium Foley Catheter		3-Way Premium Foley Catheter		Size (Ch/Fr)	Premium Foley Catheter Gold, Siliconised
3-5ml/cc	5-15ml/cc	30-45ml/cc	5-15ml/cc	30-45ml/cc			
GD23120605	-	-	-	-	-	06	Note: HV = Hard Valve
GD23120805	-	-	-	-	-	08	
GD23121005	-	-	-	-	-	10	
-	GD21111210	GD21111230	-	-	-	12	
-	GD21111410	GD21111430	-	-	-	14	
-	GD21111610	GD21111630	GD25121610	GD25121630	16		
-	GD21111810	GD21111830	GD25121810	GD25121830	18		
-	GD21112010	GD21112030	GD25122010	GD25122030	20		
-	GD21112210	GD21112230	GD25122210	GD25122230	22		
-	GD21112410	GD21112430	GD25122410	GD25122430	24		
-	GD21112610	GD21112630	GD25122610	GD25122630	26		

2-Way				Size (Ch/Fr)	2-Way Premium Hematuria Reinforced Foley Catheter Gold, Siliconised
Couvellaire		Dufour			
30-50 ml/cc	60-80 ml/cc	30-50 ml/cc	60-80 ml/cc		
GD21421830	-	-	-	18	
GD21422030	GD21422060	GD21522030	GD21522060	20	
GD21422230	GD21422260	GD21522230	GD21522260	22	
GD21422430	GD21422460	GD21522430	GD21522460	24	
GD21422630	GD21422660	GD21522630	GD21522660	26	

3-Way				Size (Ch/Fr)	3-Way Premium Hematuria Reinforced Foley Catheter Gold, Siliconised
Couvellaire		Dufour			
30-50 ml/cc	60-80 ml/cc	30-50 ml/cc	60-80 ml/cc		
GD25421830	-	-	-	18	
GD25422030	GD25422060	GD25522030	GD25522060	20	
GD25422230	GD25422260	GD25522230	GD25522260	22	
GD25422430	GD25422460	GD25522430	GD25522460	24	
GD25422630	GD25422660	GD25522630	GD25522660	26	

Catalogue No.	Size (Ch/Fr)	Nelaton Premium Foley Catheter Gold, Siliconised
GD36101200	12	
GD36101400	14	
GD36101600	16	
GD36101800	18	
GD36102000	20	
GD36102200	22	
GD36102400	24	
GD36102600	26	



## Adult Urine Collecting Bag

- Intended for even filling and complete drainage of urine.
- Eyelets allow easy hanging and provide maximum protection against tearing.
- Easy to read scale for quick determination of urine volume.
- With non-return valve.

CATALOGUE NO.	SIZE
IC1064(A)	2000 ml
IC1064SH(A)	2000 ml (with urine bag hanger)



## Paediatric Urine Collector

- For collection of urine from paediatric.
- Suitable for both male and female infants.
- Capacity: 100ml.

CATALOGUE NO.	SIZE
IC1062	100 ml

## Male External Catheter

- Designed for urine incontinence for day and night use in male patients.
- Proximal end is designed for safe and easy connection to the urine bag.

CATALOGUE NO.	SIZE
IC1010-S	Small
IC1010-M	Medium
IC1010-L	Large
IC1010-XL	Extra Large



## Urine Collecting Bag with Measured Volume Meter

- Designed for measurement of hourly urine output.
- Clear and transparent measured volume chamber provides accurate measurement.
- Urine chamber directly attached to the urine bag with overflow facility.
- Screw type bottom outlet is provided for easy and quick emptying of urine bag.
- Universal tapered connector provided sampling port facilitates midstream urine sampling.

CATALOGUE NO.	BAG CAPACITY	CHAMBER CAPACITY
IC1065	2000 ml	250 ml



- Closed urine drainage system with graduated chamber to measure the urine output on hourly basis in ICU & CCU's.
- 500 ml transparent chamber has four fine measuring sections to measure the accurate & precise volume of the urine.
- Collection chamber is provided with hydrophobic filter to minimize the risk of urinary tract infection.
- Single handed operation of the valve for emptying all chambers simultaneously.
- Bypass system is provided to prevent the overflow.
- Two litre capacity urine collecting bag with clamp & "T" shaped bottom outlet.
- Tapered connector provided sampling port facilitates easy & aseptic mid-stream urine sampling.

CATALOGUE NO.	BAG CAPACITY	CHAMBER CAPACITY
IC1066	2000 ml	500 ml

## Urine Bag Hanger

- Used for the flexible hanging of urine bags.

CATALOGUE NO.	COLOR
IC1060	Blue






## Intravenous Catheter with Injection Valve, Quick Flashback


- Quick Flashback instantly confirms successful venipuncture.
- Adva Needle Technology enhances success in first-prick of needle.
- Instant confirmation of blood flow along catheter body increases clinician's ability to successfully access the vein.



### Standard I.V. Catheter










### Quick Flashback I.V. Catheter




#### COMMON FEATURES:

- Specially designed for continuous peripheral intravenous infusion.
- Kink resistant FEP-Teflon radio-opaque catheter has extremely smooth inner and outer surface to ensure minimum trauma and smooth flow.
- Double tapered moulded catheter tip ensures smooth and less traumatic cannulation.
- Transparent flashback chamber allows immediate detection of blood.
- Proximal flexible wings for easy fixation.

CATALOGUE NO.	SIZE	COLOR CODE	OD (mm)	ID (mm)	LENGTH (mm)	FLOW RATE (ml/min)
IC1039-14	14G x 5.1cm		2.0	1.7	51	305
IC1039-16	16G x 5.0cm		1.7	1.3	50	200
IC1039-18	18G x 4.5cm		1.3	0.9	45	95
IC1039-20	20G x 3.2cm		1.1	0.8	32	65
IC1039-22	22G x 2.5cm		0.9	0.6	25	36
IC1039-24	24G x 1.9cm		0.7	0.5	19	23
IC1039-26	26G x 1.9cm		0.6	0.45	19	17

## Intravenous Catheter with Injection Port & Wings



### Sharp 3-Facet Back-Cut Needle

- Sharp and thin walled needle with reduced bevel angle for less contact area during insertion resulting **atraumatic cannulation**.

### Fluorinated Ethylene Propylene (FEP) Catheter

- Special tapered and kink resistant **Fluorinated Ethylene Propylene (FEP)** radiopaque catheter for easy insertion with optimal flow rate.

### Flexible Wing

- Special flexible wing to **provide easy and secure fixation**.

### Injection Port

- Injection port with non-return Silicone Valve for **intermittent medication**.

### Needle Grip

- Ergonomic needle grip design for **one-handed puncture technique** and comfortable handling.

### Removable Closing Cone


- Luer lock plug at the end of catheter to **prevent blood contamination**.

- ✓ Disposable.
- ✓ Sterile.
- ✓ Easy to use.
- ✓ Latex free.
- ✓ PVC free.

	CATALOGUE NO.	SIZE	COLOR CODE	OD (mm)	ID (mm)	LENGTH (mm)	FLOW RATE (ml/min)
	IC2001-14	14G	Orange	2.0	1.7	50	305
	IC2001-16	16G	Grey	1.7	1.3	50	200
	IC2001-18	18G	Green	1.3	0.9	45	95
	IC2001-20	20G	Pink	1.1	0.8	33	65
	IC2001-22	22G	Blue	0.9	0.6	25	38
	IC2001-24	24G	Yellow	0.7	0.5	19	23
	IC2001-28	26G	Purple	0.6	0.45	19	17



## Intravenous Catheter with Injection Valve



### Sharp 3-Facet Back-Cut Needle

- Sharp and thin walled needle with reduced bevel angle for less contact area during insertion resulting **atraumatic cannulation**.

### Biocompatible Polyurethane (PUR) Catheter

- Special tapered bio-compatible **Polyurethane (PUR)** catheter reduces the chances of phlebitis and thrombo-phlebitis; kink resistant allow longer indwelling time.
- Smooth surface catheter embedded with radiopaque lines for better x-ray visualization.

### Flexible Perforated Colored Wing For Ventilation

- Special flexible perforated colored wing for **secure fixation and ventilation**.

### User Friendly Snap Flip Cap & Non-Return Silicon Valve

- Injection port with snap flip cap for **one-handed operation of easy opening and closing**.
- Compatible with both luer syringe hub.
- Needle-free injection reduces risk of needlestick injury.
- Universal color-coded on port for easy gauge size identification.

### Needle Grip

- Ergonomic needle grip design for **one-handed puncture technique** and comfortable handling.

### Transparent Flashback Chamber and Hydrophobic Filter Membrane

- Transparent flashback chamber for **easy visualization** of blood.
- Hydrophobic filter to prevent blood leakage during cannulation.


### Removable Closing Cone

- Luer lock plug at the end of catheter to **prevent blood contamination**.

- ✓ Disposable.
- ✓ Sterile.
- ✓ Easy to use.
- ✓ Latex free.
- ✓ PVC free.

	CATALOGUE NO.	SIZE	COLOR CODE	OD (mm)	ID (mm)	LENGTH (mm)	FLOW RATE (ml/min)
	IC2002-14	14G	Orange	2.0	1.7	50	305
	IC2002-16	16G	Grey	1.7	1.3	50	200
	IC2002-18	18G	Green	1.3	0.9	45	95
	IC2002-20	20G	Pink	1.1	0.8	33	65
	IC2002-22	22G	Blue	0.9	0.6	25	36
	IC2002-24	24G	Yellow	0.7	0.5	19	23
	IC2002-26	26G	Purple	0.6	0.45	19	17

## Intravenous Catheter Without Port & Without Wings



### Sharp 3-Facet Back-Cut Needle

- Sharp and thin walled needle with reduced bevel angle for less contact area during insertion resulting **atraumatic cannulation**.

### Biocompatible Polyurethane (PUR) Catheter

- Special tapered bio-compatible **Polyurethane (PUR)** catheter reduces the chances of phlebitis and thrombo-phlebitis; kink resistant allow longer indwelling time.
- Smooth surface catheter embedded with **radiopaque lines** for better x-ray visualization.

### Catheter Hub

- Ergonomic product design with finger guard facilitates easy insertion **with one-handed insertion**.
- Universal color-coded on catheter hub for easy gauge size identification.

### Transparent Flashback Chamber and Hydrophobic Filter Membrane

- Transparent flashback chamber for **easy visualization** of blood.
- Hydrophobic filter to prevent blood leakage during cannulation.
- Ribbed design on chamber secure handling during cannulation.

### Removable Closing Cone

- Luer lock plug at the end of catheter to **prevent blood contamination**.

- ✓ Disposable.
- ✓ Sterile.
- ✓ Easy to use.
- ✓ Latex free.
- ✓ PVC free.

	CATALOGUE NO.	SIZE	COLOR CODE	OD (mm)	ID (mm)	LENGTH (mm)	FLOW RATE (ml/min)
	IC2003-14	14G	Orange	2.0	1.7	45	305
	IC2003-16	16G	Grey	1.7	1.3	45	200
	IC2003-18	18G	Green	1.3	0.9	45/32	90
	IC2003-20	20G	Red	1.1	0.8	32	56
	IC2003-22	22G	Blue	0.9	0.6	25	36
	IC2003-24	24G	Yellow	0.7	0.5	19	23
	IC2003-26	26G	Purple	0.6	0.45	19	17



## Blood Giving Set with Airvent (Blood Transfusion Set)

- Transparent and flexible double cylindrical drip chamber with/without ring.
- Provided with 200 micron blood filter.
- Sharp and easy piercing air vented spike.
- Soft and kink resistant tubing to ensures uniform flow rate.
- Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates.
- Self sealing 'Y' injection port designed for intermittent medications.
- Luer lock fitting with hydrophobic priming filter.
- Tube length: 180cm (130cm + 50cm).

CATALOGUE NO.	SIZE	TYPE
IC3050	Standard	Double Chamber
IC3052	Standard	Single Chamber

**Optional:**  
Available separately  
packed needle size  
18G X 1 ½"



Single Chamber



Double Chamber



Y-Connection with Protective Shield



Needle Free Y-Connection

## Intravenous Administration Set

- Transparent and flexible drip chamber.
- Provided with disc type fluid filter.
- Sharp and easy piercing air vented spike.
- Soft and kink resistant tubing to ensures uniform flow rate.
- Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates.
- Self sealing 'Y' injection port with protective shield designed for intermittent medications (Y-Connection with Protective Shield).
- Needle free 'Y' injection port with closing cap designed for intermittent medications (Needle Free Y-Connection, Sterile).
- Approximate: 20 drops / ml
- Tube length: 180cm (150cm + 30cm).

CATALOGUE NO.	SIZE	TYPE
IC3053A-21	21G X 1 ½"	Y-Connection with Protective Shield
IC3053A-23	23G X 1"	
IC3053B	Standard	Needle Free Y-Connection

## Intravenous Administration Set

- Transparent and flexible double cylindrical drip chamber with ring for secure gripping (Double Chamber & Y-Connection with Protective Shield).
- Transparent and flexible drip chamber.
- Provided with disc type fluid filter.
- Sharp and easy piercing air vented spike.
- Soft and kink resistant tubing to ensures uniform flow rate.
- Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates.
- Self sealing latex free 'Y' injection port with protective shield designed for intermittent medications.
- Luer Lock fitted with needle size 21 x 1 ½" or 23G x 1" (Double Chamber & Y-Connection with Protective Shield).
- Approximate: 20 drops / ml.
- Tube length: 150cm (120cm + 30cm).

CATALOGUE NO.	SIZE	TYPE
IC3054-21	21G x 1 ½"	Double Chamber & Y-Connection with Protective Shield
IC3054-23	23G x 1"	Double Chamber & Y-Connection with Protective Shield
IC3055-21	21G X 1 ½"	Latex Free Y-Connection
IC3055-23	23G X 1"	Latex Free Y-Connection



Double Chamber & Y-Connection with Protective Shield



Latex Free Y-Connection

**Optional:**

- Available separately packed needle size 21G X 1 ½" & 23G X 1"



Adult, 20 Drops Fitted with 21G Needle

**Optional:**

- Available with luer-lock with rubber y-site



Paediatric, 60 Drops Fitted with 23G Needle

- Transparent and flexible drip chamber.
- Provided with disc type fluid filter.
- Sharp and easy piercing air vented spike.
- Soft and kink resistant tubing to ensures uniform flow rate.
- Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates.
- Latex flashbulb injection site designed for intermittent medications.
- Luer lock fitted with needle.
- Tube length: 150cm.

CATALOGUE NO.	SIZE	TYPE
IC3056	21G X 1 ½"	Adult, 20 Drops Fitted with 21G Needle
IC3057	23G X 1"	Paediatric, 60 Drops Fitted with 23G Needle



### Microdrip Precision Solution Administration Set

- Transparent, soft, cylindrical & calibrated measured volume chamber with bold marking scale.
- Sharp piercing spike provided with built in air inlet with bacterial barrier filter.
- Automatic shut off valve prevents air embolism in the fluid line once the measured volume chamber is empty, and get as volume indicator.
- 'Y' injection port designed for additional medication.
- Approximate: 60 drops / ml.



CATALOGUE NO.	SIZE
IC3058	110 ml
IC3059	150 ml



### Peritoneal Dialysis Transfusion Set

- Specially designed 'Y' shaped transfusion set for administering dialysis solution.
- Sharp and easy piercing air vented spike.
- Provided with two upper control clamps to facilitate alternative change of fluid bottles.
- Two lower control clamps help easy input and drainage of solution.
- Proximal end is provided with suitable adapter for easy connection to catheter unit.

CATALOGUE NO.	TYPE
IC3012-10	Standard

### Peritoneal Dialysis Catheter Set

- Stainless steel trocar is provided with triple angle sharpness to facilitate smooth penetration.
- Spare scalpel blade is supplied to provide the incision.
- Multihole, open end catheter is manufactured from medical grade polyamide.
- The trocar and the tip of the catheter is perfectly matched to facilitate trauma free entry.
- Latex free flash ball type injection port is provided on the junction unit for additional medication.
- Complete kit is placed in thermoformed tray.



CATALOGUE NO.	SIZE
IC3010-25A	Adult
IC3010-25C	Child

## Extension Tube With Male Female Luer Lock

- Small bore high pressure extension line suitable for high pressure monitoring.
- Kink resistant tubing.
- Approximate prime volume: I.D: 1.0 mm & O.D: 2.0 mm.

CATALOGUE NO.	LENGTH
IC3036-200	200 cm



## CVP Manometer

- Suitable for measure and monitor continuous or intermittent monitoring of central venous pressure during infusion.
- Manometer tube with graduation from -4 cm to +34 cm and attached to three way tap.
- 120cm long extension tube is fitted with female luer lock on one end for safe connection to central venous catheter.
- Moulded clamps with latex free rubber strings are provided for easy fixation on the I.V. stand.

CATALOGUE NO.	TYPE
IC3008	Standard

## Extension Tube with 'T' Connector

- Extension Tube with T type luer lock connector.
- Intergrated with an intermittent injection port at one end, line clamp and female luer lock connector at the other end.
- Individually packed in peelable, soft blister pack.

CATALOGUE NO.	SIZE
IC1610-01	13 cm





## Ryle's Tube

- Suitable for nasogastric introduction, for nutritional purposes or aspiration of intestinal secretion.
- Manufactured from medical grade PVC.
- Super smooth, low friction surface facilitates easy intubation & reduces the incidence of nasal discharge.
- Radio opaque line is provided throughout the length.
- Distal end is coned with radio opaque material.
- Proximal end is provided with dual type universal funnel shape connector & female luer mount with closure cap.



Dual type universal funnel.

CATALOGUE NO.	SIZE	COLOR CODE
IC4032-08	08 Fg	
IC4032-10	10 Fg	
IC4032-12	12 Fg	
IC4032-14	14 Fg	
IC4032-16	16 Fg	
IC4032-18	18 Fg	
IC4032-20	20 Fg	

## Infant Feeding Tube

- Manufactured from DEHP Free PVC.
- Suitable in neonatal and paediatric enteral nutritional feeding.
- Closed distal end with two lateral eyes.
- Tube with radio-opaque line, marked at 20 cms from the tip for accurate placement.




CATALOGUE NO.	SIZE	COLOR CODE
IC4013-05	5 Fg	
IC4013-06	6 Fg	
IC4013-08	8 Fg	
IC4013-10	10 Fg	

## Umbilical Catheter

- Designed for intermittent or continuous access to the umbilical vessels of new born or pre-mature infants.
- Smooth round tip helps easy cannulation.
- Open distal end without lateral eyes, eliminates the chance of clot formation in the blind spaces.
- Radio opaque line is provided for x-ray visualisation.

SIZE	CATALOGUE NO.	COLOR CODE
3.5 Fg	IC3002-35	
4 Fg	IC3002-04	
5 Fg	IC3002-05	
6 Fg	IC3002-06	
8 Fg	IC3002-08	



## Mucus Extractor

- For aspiration of secretion from oropharynx to ensure free respiration in newly born babies.
- Convenient to extract mucus specimen for microbiological examination.
- Separate screw top lid for safe sealing and handling of aspirate.
- Soft, frosted & kink resistant PVC tubing.

CATALOGUE NO.	SIZE
IC5018-06	6 Fg
IC5018-08	8 Fg
IC5018-10	10 Fg
IC5018-12	12 Fg
IC5018-14	14 Fg



## Surgery & Wound Drainage

### Under Water Seal Pleural Drainage System



- Available in different capacities.
- Kink resistant large bore tubing facilitates unrestricted flow.
- Drainage bottle suitable for pleural drainage in conjunction with intercostal drainage catheters.

CATALOGUE NO.	SIZE
IC5020	2000 ml
IC5031	1200 ml
IC5032	500 ml

### Connecting Tube with Clamp & 'Y' Connector

- Provided with Radio-opaque Line and satin smooth eyes for drainage.
- Kink resistant.

CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
IC5009-06	6 Fg	IC5009-14	14 Fg
IC5009-08	8 Fg	IC5009-16	16 Fg
IC5009-10	10 Fg	IC5009-18	18 Fg
IC5009-12	12 Fg		



### Perforated Redon Drain Catheter

- Provided with radio-opaque line and smooth eyes.
- Soft, frosted & kink resist PVC tubing.

CATALOGUE NO.	SIZE
IC9006-08	8 Fg
IC9006-10	10 Fg
IC9006-12	12 Fg
IC9006-14	14 Fg
IC9006-16	16 Fg
IC9006-18	18 Fg



## Close Wound Suction Unit

- Suitable for close wound drainage under negative pressure post operatively with the option to operate one or two catheters simultaneously.
- Connecting tube is kink resistant and provided with additional strength to withstand the suction.
- Easy to depress chamber to activate the suction of bellow unit, single handedly.

CATALOGUE NO. (Capacity: 50 ml)	SIZE	CATALOGUE NO. (Capacity: 400-600 ml)	SIZE
IC5006-08	8 Fg	-	8 Fg
-	10 Fg	IC5008-10	10 Fg
-	12 Fg	IC5008-12	12 Fg
-	14 Fg	IC5008-14	14 Fg
-	16 Fg	IC5008-16	16 Fg
-	18 Fg	IC5008-18	18 Fg



## Suction Connecting Tube (CE certified)


- Suitable for pre-operative removal of secretion and body fluids.
- Universal connector due to flexible connectors.
- Soft and kink resistant ribbed PVC tubing.
- Individually packed perforated double wrapped packaging.

CATALOGUE NO.	LENGTH
IC551-020	2.0 m
IC551-030	3.0 m

## Yankaur Suction Set (CE certified)

- Suitable for convenient removal of secretion, blood and debris.
- Soft and kink resistant ribbed PVC tubing.
- Individually perforated double wrapped packaging.

CATALOGUE NO.	LENGTH
IC5016A	Standard



## Yankaur Suction Handle (CE certified)

- Transparent moulded handle for pre-operative suction.
- Yankaur Handle is open end with lateral eyes to withstand negative pressure and ensure continuous and easy suction of fluid.

CATALOGUE NO.	SIZE
IC9014DS	Standard tip
IC9014DC	Crown tip

## Colostomy Bag

- Designed for closed colostomy management.
- Made of medical grade PVC.
- The hole opening can be increased as per requirement.

CATALOGUE NO.	SIZE
IC4020	140mm x 250mm



## Thoracic Trocar Catheter

- Specially designed for quick non operative pleural and chest drainage for the relief of chest conditions such as tension pneumothorax.
- Blunt yet smooth tip.
- Two lateral eyes for efficient drainage.

CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
IC5030-08	8 Fg	IC5030-20	20 Fg
IC5030-10	10 Fg	IC5030-24	24 Fg
IC5030-12	12 Fg	IC5030-28	28 Fg
IC5030-14	14 Fg	IC5030-32	32 Fg
IC5030-16	16 Fg	IC5030-36	36 Fg
IC5030-18	18 Fg	IC5030-40	40 Fg

## Chest Drainage Catheter straight & Right angle

- Most suitable for post operative drainage following cardio-thoracic surgeries.
- Distal end is open with super smooth large eyes to maximise the drainage from the chest.
- Proximal end is fitted with specially designed tapered tongue to provide excellent forceps grip and non-traumatic insertion.

### Straight

CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
IC5010S-12	12 Fr	IC5010S-28	28 Fr
IC5010S-16	16 Fr	IC5010S-30	30 Fr
IC5010S-18	18 Fr	IC5010S-32	32 Fr
IC5010S-20	20 Fr	IC5010S-36	36 Fr
IC5010S-22	22 Fr	IC5010S-40	40 Fr
IC5010S-24	24 Fr	-	-

### Right Angle

- Provided with curved (90°) angle for easy aspiration of drainage.

CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
IC5010A-12	12 Fr	IC5010A-24	24 Fr
IC5010A-16	16 Fr	IC5010A-28	28 Fr
IC5010A-20	20 Fr	IC5010A-32	32 Fr





## Petri Dishes

- Made of clear general purpose polystyrene.
- Easy to stack and grip.
- Suitable for machine operating.

CATALOGUE NO.	SIZE
IC24-3001R	90 mm x 15 mm (Single compartment)



\* Patented

Trunk End

Fluted-beaked

## Proctoscope

- Made from clear and rigid clear General Purpose polystyrene.
- Used in diagnosis of rectal examination.
- The device has two parts: obturator and speculum.
- Individually pack in polybag.

CATALOGUE NO.	SIZE	TYPE
IC8900A	Standard	Trunk End
IC8900	Standard	Fluted-beaked

## Urine Specimen Container

- Specifically designed for collection, transport, and storage of specimens.
- Containers are translucent for easy viewing contents.
- Available in various color caps.
- Also available in sterile and non-sterile.

POLYPROPYLENE		POLYSTYRENE	
CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
IC1061RPP-NS	60 ml	IC1063RPS-S	60 ml
IC1061WPP-NS	60 ml	IC1063WPS-S	60 ml
IC1061YPP-NS	60 ml	IC1063YPS-S	60 ml
IC1062SC-90	90 ml		



## Sputum Cup


- Circular flat shape.
- Screw type circular lid with smooth upper surface.

CATALOGUE NO.	SIZE
IC2054-500	2 oz

## Amniotic Membrane Perforator

- Anatomically shaped to allow user to snag membrane, preventing injury to the fetus.
- Curved, protected hook, to safely rupture the amniotic membrane.

CATALOGUE NO.	SIZE
IC2052S (Straight)	25 cm



## Vaginal Speculum

- Side screw type with various sizes.
- Strong and simple to use.
- Screw thread ratchet system with nut, patient friendly.

CATALOGUE NO.	SIZE
IC2043-S	Small
IC2043-M	Medium
IC2043-L	Large

## Cervical Brush

- Used for cervical cell samplings.

CATALOGUE NO.	LENGTH
IC6026	19 cm



## Gynecological Brush with Tip Ball End

- Used for endocervical samplings.



CATALOGUE NO.	LENGTH
IC6025-100	19 cm

## Procedure & Clinical Sets

### Basic Dressing Set (Premium)

- 1 pce Tray of 3 Compartments
- 2 pcs Forceps (Blue)
- 1 pce Forcep (Red)
- 5 pcs Cotton Balls
- 5 pcs Gauze Swabs (7.5 cm x 7.5 cm - 8ply)
- 1 pce C-Fold Hand Towel
- 1 pce PE Patient Drape (45 cm x 45 cm)
- 1 pce Yellow Limpet Bag
- 1 pce Sterile Field-Outer Wrap

CATALOGUE NO.	TYPE
IC7000P108	3 Forceps





### Basic Dressing Set (Standard)

- 1 pce Tray of 3 Compartments
- 2 pcs Forceps
- 5 pcs Cotton Balls
- 4 pcs Gauze Swabs (7.5 cm x 7.5 cm - 8ply)
- 1 pce C-Fold Hand Towel
- 1 pce PE Patient Drape
- 1 pce Yellow Limpet Bag
- 1 pce Sterile Field-Outer Wrap



CATALOGUE NO.	TYPE
IC7000S	2 Forceps



### Basic Dressing Set (Haemo)

- 1 pce Tray of 3 Compartments
- 1 pce Forcep
- 6 pcs Cotton Balls
- 4 pcs Gauze Swabs (7.5 cm x 7.5 cm - 8ply)
- 1 pce Yellow Limpet Bag
- 1 pce C-Fold Hand Towel
- 1 pce PE Patient Drape
- 1 pce Sterile Field-Outer Wrap

CATALOGUE NO.	TYPE
IC7000H	1 Forcep

### Vaginal Examination Set

- 1 pce Gallipot (250 ml)
- 1 pce Yellow Limpet Bag
- 1 pce C-Fold Hand Towel
- 2 pcs Gauze Swabs (7.5 cm x 7.5 cm - 8ply)
- 1 pair Latex Examination Glove
- 10 pcs Cotton Balls
- 1 pce Sanitary Pad with Loop
- 1 pce Sterile Field-Outer Wrap
- 1 pce Patient Drape

CATALOGUE NO.	TYPE
IC7001D	Loop
IC7001D-1	adhesive pad



### Catheterization Set

- 1 pce Kidney Tray
- 1 pce Tray
- 1 pce Sterile Field-Outer Wrap
- 1 pce Non Woven Drape
- 1 pair Latex Examination Glove
- 1 pce Fenestrated Drape (46.6 cm x 60 cm)
- 6 pcs Cotton Balls
- 2 pce Gauze Swabs (7.5 cm x 7.5 cm - 8ply)
- 2 pcs Forceps (Green)
- 1 pce Syringe 30 ml w/o Needle
- 2 pce C-Fold Hand Towel
- 1 pce Specimen Container 60ml With Label
- 1 pce Yellow Limpet Bag
- 1 pce Lubricating Jelly

CATALOGUE NO.	TYPE
IC7002D	Standard

## Delivery Set

- 1 pce Sterile Field-Outer Wrap
- 1 pair Latex Examination Glove
- 1 pce C-Fold Hand Towel
- 2 pcs Clinical Waste Bags
- 1 pce Kidney Tray
- 1 pce Tray
- 1 pce Under pad / Liner (Patient Drape 18"x24")
- 3 pcs Sanitary Pads
- 1 pce Umbilical Cord Clamp
- 30 pcs Cotton Balls
- 1 pce Tape
- 1 pce Tampon

CATALOGUE NO.	TYPE
IC7003D	Standard



## Pap Smear Test Kit

- 1 pce Gynecological Brush with Tip Ball End
- 1 pce Cervical Brush
- 1 pce Glass Slide with Frosted End
- 1 pce Plastic Glass Slide Holder
- 1 pce Patient Indication Label

CATALOGUE NO.	TYPE
IC7006D	Standard



## Top & Tail / Baby Cleaning Set

- 1 pce Gallipot (450 ml)
- 15 pcs Cotton Balls
- 1 pce Sterile Field-Outer Wrap

CATALOGUE NO.	TYPE
IC7008D	Standard



## Injection Set

- 1 pce Tray of 3 compartments
- 2 pcs Forceps
- 1 pce Gauze Swabs (7.5 cm x 7.5 cm - 8ply)
- 3 pcs Cotton Balls
- 1 pce Yellow Limpet Bag
- 1 pce Sterile Field-Outer Wrap
- 1 pce C-Fold Hand Towel

CATALOGUE NO.	TYPE
IC7005D	Standard

## Toilet & Suture Set

- 1 pce Gallipot Tray 950ml, Rectangle with Cover
- 2 pcs Cotton Balls
- 2 pcs Forceps
- 6 pcs Gauze Swabs (7.5 cm x 7.5 cm - 8ply)
- 1 pce Underpad (17" x 24")
- 1 pce Crepe Paper (500 mm x 500 mm)
- 1 pair Latex Examination Glove (Medium)
- 2 pcs C-Fold Hand Towels
- 1 pce Yellow Limpet Bag

CATALOGUE NO.	TYPE
IC7007D	Standard



### Eye Dressing Set

- 1pce Tray of 3 compartments
- 3pcs Forceps
- 6pcs Cotton Balls
- 4pcs Gauze Swabs (7.5 cm x 7.5 cm - 8ply)
- 1pce C-Fold Hand Towel
- 1pce Plastic Patient Drape
- 1pce Yellow Limpet Drape
- 1pce Sterile Field-Outer Wrap
- 2pcs Adhesive Eye Pads



CATALOGUE NO.	TYPE
IC7010D	Standard



### Cord Dressing Set

- 1pce Plastic Gallipot (60 ml)
- 5pcs Cotton Wool Balls (0.6 gm)
- 1pce Forcep (Blue)
- 1pce Yellow Limpet Bag

CATALOGUE NO.	TYPE
IC7011D-3	Standard



### Suction Set

- 1 pce Gallipot (120 ml)
- 3 pcs Gauze Swabs (7.5 cm x 7.5 cm - 8ply)
- 1 pce Sterile Field-Outer Wrap

CATALOGUE NO.	SIZE
IC7013D	Standard



### Blood Taking Set

- 1 pce Tray of 3 compartments
- 1 pce Adhesive Plaster
- 1 pce PE Patient Drape
- 1 pce Sterile Field-Outer Wrap
- 2 pcs Gauze Swabs (5 cm x 5 cm - 8ply)
- 4 pcs Cotton Balls

CATALOGUE NO.	TYPE
IC7020	Standard

### Swabbing Set

- 1 pce Gallipot (450 ml)
- 20 pcs Cotton Balls

CATALOGUE NO.	TYPE
IC7021B	Standard





## Kidney Dish

- Used in medical and surgical wards to receive soiled dressings and other medical waste.
- Optional: Available in various color (yellow, red, green, blue and white)

CATALOGUE NO.	SIZE
IC2020-PB	700 ml

## Forceps

- For grasping or holding objects.
- Control mechanism for pressure application.
- Optional: Available in various color.

CATALOGUE NO.	COLOR
IC6029B	Blue



## Umbilical Cord Clamp

- Suitable for clamping the umbilical cord of new born baby immediately after birth.

CATALOGUE NO.	SIZE
IC5084	Standard

# Wound Care & Dressing

## Cotton Wool Roll

CATALOGUE NO.	SIZE
IC2158-500	500 grams



## Combine Dressing Roll

CATALOGUE NO.	SIZE
IC2159-500	500 grams (20 cm x 10 M)

## White Elastic Adhesive Bandage

CATALOGUE NO.	SIZE
IC640275-075	7.5 cm x 2.75 M





**Gauze Ribbon**

CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
IC2093-122024	1.25 cm x 5.5 M	IC2093-902024	5 cm x 90 M
IC2093-551915	5 cm x 5.5 M	IC2093-102024	2.5 cm x 100 M



**Lint Absorbent**

CATALOGUE NO.	SIZE
IC2160-500	500 grams

**Plaster Extension**

CATALOGUE NO.	SIZE
IC093003-02	10 cm x 2.5 M



**Skin Color Crepe Bandage**

CATALOGUE NO.	SIZE
IC2121-050	5 cm x 4 M
IC2121-075	7.5 cm x 4 M
IC2121-100	10.0 cm x 4 M
IC2121-150	15.0 cm x 4 M

**Cotton Wool Balls**

CATALOGUE NO.	SIZE
IC2151-005A	0.5 grams



## Transparent Film Dressing

CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
IC2162-010	10 cm x 12 cm	IC2162-065	6.5 cm x 5 cm
IC2162-010B	10 cm x 20 cm	IC2162-095	9.5 cm x 8.5 cm
IC2162-010C	10 cm x 25 cm	IC2162-155	15.5 cm x 8.5 cm
IC2162-040	4 cm x 5 cm	IC2162-235	23.5 cm x 8.5 cm
IC2162-060	6 cm x 7 cm		



## Paraffin Gauze Dressing

CATALOGUE NO.	SIZE
IC5.02	10 cm x 10 cm



## Transparent Infusion Set Dressing

CATALOGUE NO.	SIZE
IC2163-067	6 cm x 7 cm



## Soft Paper Surgical Tape

CATALOGUE NO.	SIZE
IC2110-500	5.0 cm x 9.1 M
IC2110-750	7.5 cm x 9.1 M

## Surgical Drape

CATALOGUE NO.	SIZE
IC2048-270	270 cm x 93 cm
IC2048-38	38 cm x 45 cm







### Alcohol Swabs

CATALOGUE NO.	SIZE
IC0645	2 ply-30 x 30mm (folded)

### Wooden Applicator (Orange Stick)

CATALOGUE NO.	SIZE
IC2155	15 cm



### Wooden Tongue Depressor

CATALOGUE NO.	SIZE
IC2154	150 mm x 18 mm x 2 mm

### Wheelchair

CAT. NO.	SIZE
IC869X-46	Standard

**Available in:**  
Solid & Pneumatic Tyre




## Rehabilitation





## Our Accreditation



**DNV GL Certification  
Body, Norway**



**Kgs Certification Sdn. Bhd.  
Malaysia**



**Ministry Of Health,  
Malaysia**



**Medical Device Authority,  
Malaysia**

## Our Achievement



**Enterprise 50 Award  
(champion)**



**FMM Icon Award,  
SME Category**



**AKI Award,  
SME Category**



**MATRADE Export Award,  
Bumiputera Export Category**

 [www.idealcare.com.my](http://www.idealcare.com.my)



**IDEALCARE®**  
.....the quality you can trust

**Head Office / Distribution Centre:**

No. 70 & 71, Jalan Sungai Tukang 2/1,  
Kawasan Perusahaan Sungai Tukang,  
08000 Sungai Petani, Kedah, Malaysia.  
Tel : 6 (04) - 422 7775, (6) 04 - 422 8775  
Fax : 6 (04) - 422 6775  
E-mail : [info@idealcare.com.my](mailto:info@idealcare.com.my)

**Center of Excellence / COE**

Lot PT. 60598,  
Kawasan Perindustrian LPK Fasa 1,  
08000 Sungai Petani, Kedah, Malaysia.  
Tel : 6 (04) - 444 6875, (6) 04 - 444 6975  
Fax : 6 (04) - 444 5975  
E-mail : [info@idealcare.com.my](mailto:info@idealcare.com.my)