

**INFUSION / TRANSFUSION  
THERAPY**

**IC3050 (17.0601.01)  
Blood Giving Set with Airvent  
(Blood Transfusion / Intravenous Admin Set,  
Double Chamber), Sterile**

- Transparent and flexible double cylindrical drip chamber with ring.
- Provided with 200 micron blood filter.
- Sharp and easy piercing air vented spike.
- Soft and kink resistant tubing to ensures uniform flow rate.
- Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates.
- Self sealing 'Y' injection port designed for intermittent medications.
- Luer lock fitting with hydrophobic priming filter.
- Tube length: 180cm (130cm + 50cm).



**IC3052  
Blood Giving Set with Airvent  
(Blood Transfusion / Intravenous Admin Set), Sterile**

- Transparent and flexible cylindrical drip chamber.
- Provided with 200 micron blood filter.
- Sharp and easy piercing air vented spike.
- Soft and kink resistant tubing to ensures uniform flow rate.
- Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates.
- Self sealing 'Y' injection port designed for intermittent medications.
- Luer lock fitting with hydrophobic priming filter.
- Tube length: 180cm (130cm + 50cm).



- Optional:**  
• Available separately packed needle size 18G X 1 1/2"



**IC3053A (17.0602.01)  
Intravenous Admin Set with Airvent & Y-Connection  
with Protective Shield, Sterile**

- Transparent and flexible drip chamber.
- Provided with disc type fluid filter.
- Sharp and easy piercing air vented spike.
- Soft and kink resistant tubing to ensures uniform flow rate.
- Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates.
- Self sealing 'Y' injection port with protective shield designed for intermittent medications.
- Approximate: 20 drops / ml
- Tube length: 180cm (150cm + 30cm).



- Optional:**  
• Available separately packed needle size 21G X 1 1/2" & 23G X 1"



**IC3053B  
Intravenous Admin Set with Airvent & Needle Free  
Y-Connection, Sterile**

- Transparent and flexible drip chamber.
- Provided with disc type fluid filter.
- Sharp and easy piercing air vented spike.
- Soft and kink resistant tubing to ensures uniform flow rate.
- Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates.
- Needle free 'Y' injection port with closing cap designed for intermittent medications.
- Approximate: 20 drops / ml.
- Tube length: 180cm (150cm + 30cm).



**IC3054  
Intravenous Admin Set Double Chamber with Airvent &  
Y-Connection with Protective Shield, Sterile**

- Transparent and flexible double cylindrical drip chamber with ring for secure gripping.
- Provided with disc type fluid filter.
- Sharp and easy piercing air vented spike.
- Soft and kink resistant tubing to ensures uniform flow rate.
- Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates.
- Self sealing latex free 'Y' injection port with protective shield designed for intermittent medications.
- Luer Lock fitted with needle size 21 x 1 1/2" or 23G x 1".
- Approximate: 20 drops / ml.
- Tube length: 150cm (120cm + 30cm).



**IC3055  
Intravenous Admin Set with Airvent & Latex Free  
Y-Connection, Sterile**

- Transparent and flexible drip chamber.
- Provided with disc type fluid filter.
- Sharp and easy piercing air vented spike.
- Soft and kink resistant tubing to ensures uniform flow rate.
- Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates.
- Self sealing latex free 'Y' injection port designed for intermittent medications.
- Approximate: 20 drops / ml
- Tube length: 150cm (120cm + 30cm).



- Optional:**  
• Available separately packed needle size 21G X 1 1/2" & 23G X 1"



**IC3056  
Intravenous Admin Set with Airvent Adult, 20 Drops  
Fitted with 21G Needle, Sterile**

- Transparent and flexible drip chamber.
- Provided with disc type fluid filter.
- Sharp and easy piercing air vented spike.
- Soft and kink resistant tubing to ensures uniform flow rate.
- Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates.
- Latex flashbulb injection site designed for intermittent medications.
- Luer lock fitted with needle size 21G x 1 1/2".
- Approximate: 20 drops / ml.
- Tube length: 150cm.



**IC3057  
Intravenous Admin Set with Airvent Paediatric, 60 Drops  
Fitted with 23G Needle, Sterile**

- Flexible, soft and transparent drip.
- Provided with disc type fluid filter.
- Sharp and easy piercing air vented spike.
- Soft and kink resistant tubing to ensures uniform flow rate.
- Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates.
- Latex flashbulb injection site designed for intermittent medications.
- Luer lock fitted with needle size 23G x 1".
- Approximate: 60 drops / ml.
- Tube length: 150cm.





**DEHP  
Free**

**IC3058 (110ml) (17.0604.01)**

**IC3059 (150ml)**

**Microdrip Precision Solution Administration Set, Sterile**

- Transparent, soft, cylindrical & calibrated measured volume chamber with bold marking scale.
- Sharp piercing spike provided with built in air inlet with bacterial barrier filter.
- Automatic shut off valve prevents air trapping in the fluid line once the measured volume chamber is empty, and it also indicates the level of liquid.
- Soft and kink resistant tubing.
- 15 µm fluid filter to reduce particulate matter.
- 'Y' injection port designed for additional injection.
- Approximate: 60 drops / ml.



**IC3012-10 (17.0004.01)**

**Peritoneal Dialysis Transfusion Set**

**(Administration Set for Peritoneal Dialysis), Sterile**

- Specially designed 'Y' shaped transfusion set for administering dialysis solution.
- Sharp and easy piercing air vented spike.
- Provided with two upper control clamps to facilitate alternative change of fluid bottles.
- Two lower control clamps help easy input and drainage of solution.
- Proximal end is provided with suitable adapter for easy connection to catheter unit.



**IC3010-25 (17.3802.01)**

**Peritoneal Dialysis Catheter Set**

**(Trocar Cannula for P.D.S.), Sterile**

- Stainless steel trocar is provided with triple angle sharpness to facilitate smooth penetration.
- Spare scalpel blade is supplied to provide the incision.
- Multihole, open end catheter is manufactured from medical grade polyamide.
- The trocar and the tip of the catheter is perfectly matched to facilitate trauma free entry.
- Latex free flash ball type injection port is provided on the junction unit for additional medication.
- Complete kit is placed in thermoformed tray.

**Optional:**

- Available size: Adult/Child