

# DISPOSABLE PROTECTIVE COVERALL

PERSONAL PROTECTIVE EQUIPMENT

**IDEALCARE®**

.....the quality you can trust

**Sealing with rubber strip**

- Effective isolation suture area zipper and other locations bacterial virus invasion.

**Smooth surface without moisture**

- Can effectively prevent bacteria from adhering.

**Hood cover.**

**Waist elastic design**

- To meet the needs of different figures.

**Elastic mouth design**

- The cuff trouser leg is designed with elastic opening for easy expansion.

**Features**

- Preventing infiltration of airbornes particles and body fluid while attending to patients in healthcare facility.
- Comfortable and breathable.
- Full body coverage.

**Material**

- Polypropylene (PP) with Polyethylene (PE).
- Weight: 65 gsm.

**NON-STERILE**

**INDIVIDUAL PACKED**

Catalogue No.	Size
IC2108-165	M
IC2108-170	L
IC2108-175	XL
IC2108-180	XXL
IC2108-185	XXXL

MDR-20210125-26872

"Your One Stop Centre for Medical & Healthcare Products"



Ideal Health care Sdn. Bhd. (484861-V)  
 No. 70 & 71, Jalan Sungai Tukang 2/1, Kawasan Perusahaan Sungai Tukang,  
 08000 Sungai Petani, Kedah, Malaysia.  
 Tel : 6 (04) – 422 7775 / 422 8775 Fax : 6 (04) – 422 6775  
 Email : info@idealcare.com.my

E. L. NO.: MDA-1080-KD120

[www.idealcare.com.my](http://www.idealcare.com.my)



IHSB-CAT-BDD-PPE-11-00