

# PROTECTIVE COVERALL SUIT

PERSONAL PROTECTIVE EQUIPMENT

**IDEALCARE®**

.....the quality you can trust



**Sealing with rubber strip**

- Effective isolation suture area zipper and other locations bacterial virus invasion.

**Sealing with rubber strip**

- Effective isolation suture area zipper and other locations bacterial virus invasion.

**Features**

- Preventing infiltration of airbornes particles and body fluid while attending to patients in healthcare facility.
- Comfortable and breathable.
- Full body coverage.

**Material**

- Polypropylene (PP) with Polyethylene (PE).
- Weight: 58 gsm.
- Individual.
- Non-sterile.

**\*STERILE**

**INDIVIDUAL PACKED**

• Boot cover.




• Hooded jacket



**Waist elastic design**

- To meet the needs of different figures.



**Elastic mouth design**

- The cuff trouser leg is designed with elastic opening for easy expansion.

Catalogue No.	Size
IC2104-L (Protective coverall suit)	170 (L)
IC2104-XXL (Protective coverall suit)	180 (XXL)
IC2106 (Boot cover)	Standard

MDR-20210125-26872

"Your One Stop Centre for Medical & Healthcare Products"



Ideal Health care Sdn. Bhd. (484861-V)  
 No. 70 & 71, Jalan Sungai Tukang 2/1, Kawasan Perusahaan Sungai Tukang,  
 08000 Sungai Petani, Kedah, Malaysia.  
 Tel : 6 (04) – 422 7775 / 422 8775 Fax : 6 (04) – 422 6775  
 Email : info@idealcare.com.my

E. L. NO.: MDA-1080-KD120

[www.idealcare.com.my](http://www.idealcare.com.my)



IHSB-CAT-BDD-PPE-12-00