



3 PLY SURGICAL FACE MASK

PPE COVID 19

Components

1. Nose bar.
2. Tie-on-type & Earloop-type.
3. Ultrasonic seam.
4. Layers of medical gauze & bacterial filter.
5. Imprinted **IDEALCARE**® logo.



Features

- BFE ≥ 95% for effective particles and bacterial isolation.
- Perfect fitting, odourless, clean and highly breathable.
- Covering nose and mouth area to prevent bacteria and dirt in pandemic or toxic environment.

Manufacturer Information

IDEAL HEALTHCARE SDN. BHD.

- Establishment License No: MDA-1080-KD120.
- MDA Registration No.: GMD96851528118A.
- Shelf Life: 5 years from the manufacturing date.

New

***Also available in 4ply.**







Quality Certification



Green

Product Specification

Type	TIE-ON	EARLOOP
Product code	IC2102-003TG	IC2102-003EG
Illustration for application		
Packaging	50 pcs / box X 50 boxes /carton	
Packaging illustration		



Ideal Healthcare Sdn. Bhd.(484861-V)
No. 70 & 71, Jalan Sungai Tukang 2/1, Kawasan Perusahaan Sungai Tukang,
08000 Sungai Petani, Kedah, Malaysia.
Tel : 6 (04) – 422 7775 / 422 8775 Fax : 6 (04) – 422 6775
Email : info@idealcare.com.my

E. L. NO.: MDA-1080-KD120

"Your One Stop Centre for Medical & Healthcare Products"

www.idealcare.com.my

COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 13485 =



idealcareMY

IHSB-CAT-BDD-PPE-02-00