

TRACHEOSTOMY TUBE

ANAESTHESIA

Cuffed

Plain Uncuff

Smoothly connected inflation tube maximizes patient comfort.

Tube size and batch number are clearly marked on the pilot balloon for easy identification and traceability.

Smooth stylet minimizes tissue drag during intubation.

Soft, thin-walled and optimum sized HVLP or SVLP cuffs for effective sealing of the trachea.

Swivel 15mm connector provides easy handling.

Clear size marking on flexible neckstrap enables easy size identification and provides maximum patient comfort.

Clear radiopaque line facilitates verification of intubated position.

Preconnected necktape eliminates extra step especially during emergency situation.

Being light in weight and non irritant, it avoids inflammation around tracheostomy tube.



**Tracheostomy Tube Cuffed,
Sterile**

CATALOGUE NO.	SIZE
IC62520030	3.0 mm
IC62520035	3.5 mm
IC62520040	4.0 mm
IC62520045	4.5 mm
IC62520050	5.0 mm
IC62520055	5.5 mm
IC62520060	6.0 mm
IC62520065	6.5 mm
IC62520070	7.0 mm
IC62520075	7.5 mm
IC62520080	8.0 mm
IC62520085	8.5 mm
IC62520090	9.0 mm
IC62520095	9.5 mm
IC62520100	10.0 mm
Packaging: 10 pcs / box, sterile	

**Tracheostomy Tube Plain
Uncuff, Sterile**

CATALOGUE NO.	SIZE
IC62540030	3.0 mm
IC62540035	3.5 mm
IC62540040	4.0 mm
IC62540045	4.5 mm
IC62540050	5.0 mm
IC62540055	5.5 mm
IC62540060	6.0 mm
IC62540065	6.5 mm
IC62540070	7.0 mm
IC62540075	7.5 mm
IC62540080	8.0 mm
IC62540085	8.5 mm
IC62540090	9.0 mm
IC62540095	9.5 mm
IC62540100	10.0 mm
Packaging: 10 pcs / box, sterile	